

STERILIZATION LOG – PERSONAL SERVICE SETTINGS (PSS)

This form is for documentation of sterilization of re-usable equipment and instruments in Personal Service Settings as required in the Personal Service Settings Regulation O. Reg 136/18. For information on sterilization requirements, refer to Public Health Ontario’s Guide to Infection Prevention and Control in Personal Service Settings or contact the Health Unit.

Name of Sterilizer:						Type of Sterilization:						
Date	Time	Load No.	Physical monitoring			Chemical Monitoring (Pass/Fail)			Biological Monitoring spore test (Pass/Fail)	Instrument / Equipment	Maintenance / Repairs	Initials
			Temperature	Pressure	Length of sterilization cycle	External Indicator	Internal Indicator	Bowie-Dick Test				
						Type 1 <input type="checkbox"/> 2 <input type="checkbox"/>	Type 4 <input type="checkbox"/> 5 <input type="checkbox"/>	Yes <input type="checkbox"/> <input type="checkbox"/> n/a				